

# A NEW FRONTIER IN REGIONAL ANESTHESIA

COMBINED SPINAL-EPIDURAL KITS AND SYSTEMS

A NEW FRONTIER IN ANESTHESIA PRACTICE

## ELDOR™ Needle Minipack


Eldor Needle combined Spinal-Epidural Anesthesia



Double-Hole pencil point Spinal Needle.



Combined end - multiple lateral holes Epidural Catheter.



## ELDOR™ NEEDLE

*Combined Spinal-Epidural Needle*



- ★ A specialized needle for the combined spinal-epidural anaesthesia.
- ★ Epidural catheter can be inserted before the spinal anaesthetic injection.
- ★ Facilitates the insertion of very small gauge spinal needles through its spinal conduit, so significantly reduces the incidence of post-dural puncture headache.
- ★ No danger of epidural catheter protrusion through the dural hole made by the spinal needle.
- ★ No metallic particles production while the spinal needle passes through the bent epidural needle tip, as in the needle-through-needle technique.
- ★ After you use it successfully you will not want to use any other needle or any other technique.

**CSEN International Ltd.**

# ELDOR COMBINED-SPINAL EPIDURAL

**CSEN 7:** - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal Needle 26G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”, packed, sterilized, CE marked.

**CSEN 7a:** - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 90 mm, Eldor Spinal Needle 26G, 90 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”, packed, sterilized, CE marked.

**CSEN 7A:** - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal Needle 26G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Assist Guide, LOR Syringe. Packed, sterilized, CE marked.



# ELDOR COMBINED-SPINAL EPIDURAL

**CSEN 7B:** - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, CSEN Spinal Needle 25G, 123 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Assist Guide, LOR Syringe. Packed, sterilized, CE marked.

**CSEN 7C:** - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm; Eldor Spinal Needle 25G, 125 mm; Eldor Epidural Catheter; Adapter; Filter; Thread Assist Guide; LOR Syringe. Packed, sterilized, CE marked.

**CSEN 7D:** - Eldor Combined Spinal Epidural Kit: - Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit 120 mm; -CSEN Pencil Point Spinal Needle , Sharp tip, 25G, 125 mm, - Eldor Epidural Catheter 20G, 6 lateral holes, closed end, - Tuohy-Borst Adapter, - Epidural Catheter Flat Filter, 0.2  $\mu\text{m}$ , - Thread Assist Guide, - LOR (Loss of Resistance) Syringe 6ml, packed, sterilized, CE marked.



# ELDOR COMBINED-SPINAL EPIDURAL

**CSEN 15:** - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal Needle 27G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”, packed, sterilized, CE marked.

**CSEN 21:** - Eldor/Tuohy Combined Spinal Epidural Kit: Tuohy Needle 16G, 90 mm, Eldor Spinal Needle 26G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”, packed, sterilized, CE marked.

**CSEN 56:** - Combined Spinal Epidural Kit: Tuohy Needle 18G, 90 mm, Spinal Needle 26G, 125 mm, Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked



# ELDOR COMBINED-SPINAL EPIDURAL

**CSEN 69:** - Eldor/Tuohy Combined Spinal Epidural Kit: - Tuohy Needle 18G, 90 mm, - Eldor Spinal Needle 26G, 125 mm,- Eldor Epidural Catheter, -Adapter, - Filter, - Thread Guide, - LOR Syringe, packed, sterilized, CE marked.

**CSEN 69a:** - Eldor/Tuohy Combined Spinal Epidural Kit: - Tuohy Epidural Needle 18G, 80 mm, metal stylet, - CSEN Pencil Point Spinal Needle , Sharp tip, 25G, 123 mm, - Eldor Epidural Catheter 20G, 6 lateral holes, closed end, -Tuohy-Borst Adapter, - Epidural Catheter Flat Filter, 0.2  $\mu\text{m}$ , - CSEN Epidural Catheter Fixation Device, - Thread Assist Guide,- LOR (Loss of Resistance) Syringe 10 ml, packed, sterilized, CE marked.

**CSEN 69b:** - Eldor/Tuohy Combined Spinal Epidural Kit: - Tuohy Epidural Needle 18G, 80 mm, metal stylet, - CSEN Pencil Point Spinal Needle 25G, 123 mm, - Eldor Epidural Catheter 20G, 6 lateral holes, closed end, - Tuohy-Borst Adapter, - Epidural Catheter Flat Filter, 0.2  $\mu\text{m}$ , - Thread Assist Guide, -LOR (Loss of Resistance) Syringe 6 ml. Packed, sterilized, CE marked.

**CSEN 69c:** - Eldor/Tuohy Combined Spinal Epidural Kit: - Tuohy Epidural Needle 18G, 80 mm, metal stylet, - CSEN Pencil Point Spinal Needle , Sharp tip, 25G, 123 mm, - Eldor Epidural Catheter 20G, 6 lateral holes, closed end, -Tuohy-Borst Adapter, - Epidural Catheter Flat Filter, 0.2  $\mu\text{m}$ , - Thread Assist Guide, - LOR (Loss of Resistance) Syringe 6 ml, packed, sterilized, CE marked.



# ELDOR COMBINED-SPINAL EPIDURAL

**CSEN 69.1:** - Eldor/Tuohy Combined Spinal Epidural Kit: - Tuohy Needle 18G, 80 mm, - Eldor Spinal Needle 25G, 123 mm, - Eldor Epidural Catheter, -Adapter, - Filter, - Thread Guide, - LOR Syringe, - Fenestrated Cover 45 x70 cm with hole 10 cm, edge adhesive tape; - Syringes 2.5 ml/ 10 ml/ 20 ml; -Swab (3x); Clamp; - Compress 10 x 10 cm, 8 layers (3x); - Cannulas OD 0.5mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm;- Tray for the products; - Cover 70 x 90 cm; packed, sterilized, CE marked.

**CSEN 69.2:** CSEN Combined Spinal Epidural Kit: - Tuohy Epidural Needle 18G, 80 mm, metal stylet, - CSEN Pencil Point Spinal Needle 27G, 125 mm, -Epidural Catheter 20G, 3 lateral holes, closed end, - Tuohy-Borst Adapter, -Epidural Catheter Flat Filter, 0.2 µm, - Thread Assist Guide, - LOR Syringe 6ml. Packed, sterilized, CE marked.

**CSEN 69.3:** - Eldor/Tuohy Combined Spinal Epidural Kit: - Tuohy Needle 18G, 80 mm, - Eldor Pencil Point Spinal Needle 25G, 125 mm, - Eldor six holes Epidural Catheter, - Adapter, - Filter, - Thread Guide, - LOR Syringe, packed, sterilized, CE marked.

**CSEN 69.4:** - Eldor/Tuohy Combined Spinal Epidural Kit: - Tuohy Needle 18G, 80 mm, - CSEN Pencil Point Spinal Needle 25G, 125 mm, - CSEN three holes Epidural Catheter, - Adapter, - Filter, - Thread Guide, - LOR Syringe, packed, sterilized, CE marked.



# ELDOR COMBINED-SPINAL EPIDURAL

**CSEN 70:** - Eldor Combined Needles: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal Needle 26G, 125 mm, packed, sterilized, CE marked.

**CSEN 71:** - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal Needle 27G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

**CSEN 73:** - Eldor/Tuohy Combined Spinal Epidural Kit: Tuohy Needle 16G, 90 mm, Eldor Spinal Needle 26G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.





# ELDOR COMBINED-SPINAL EPIDURAL

**CSEN 120:** - CSEN Spinal/Epidural Tray - Polycarbonate LOR syringe 10 ml; Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x); Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

**CSEN 121:** - CSEN Spinal/Epidural Tray - LOR syringe 10 ml; Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x); Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

**CSEN 122:** - CSEN Spinal/Epidural Tray - Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x); Cannulas OD 0.5mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.



# ELDOR COMBINED-SPINAL EPIDURAL

**CSEN 127:** - Eldor Combined Spinal Epidural Tray: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal Needle 26G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”; Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x); Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

**CSEN 128:** - Eldor Combined Spinal Epidural Tray: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal Needle 26G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe; Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x); Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90cm; packed, sterilized, CE marked.

**CSEN 202:** - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal Needle 25G, 123 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

**CSEN 203:** - Eldor/Tuohy Combined Spinal Epidural Kit: Tuohy Epidural Needle 18G, 80 mm, metal stylet, CSEN Pencil Point Spinal Needle , Sharp tip, 25G, 123 mm, Eldor Epidural Catheter 20G, 6 lateral holes, closed end, Tuohy-Borst Adapter, Epidural Catheter Flat Filter, 0.2 µm, CSEN Epidural Catheter Fixation Device, Thread Assist Guide, LOR (Loss of Resistance) Syringe 10 ml, packed, sterilized, CE marked.



**THANK YOU FOR YOUR BUSINESS**

**ELDOR CSEN**

The Global Regional Anesthesia Distribution

**Distributed by**



**Focused Medical Technologies, LLC**

Focused Medical Technologies, LLC