



The Global Regional Anesthesia Distribution

CSEN International Ltd., P.O.Box 27476, Jerusalem 9127302, Israel, Tel: 972-2-5835528,
Fax: 972-2-5835744 , <http://www.csen.com/anesthesia>, <http://www.csen.com> ,
csen_international@csen.com

ELDOR CSEN Catalogue: Wound Catheters

CSEN 194: - **ELDOR Wound Catheter kit** - OD 1 mm x L 1000 mm, 10

lateral holes, 10 mm apart, closed tip; Wound catheter adapter (Tuohy Borst);

Wound catheter thread assist guide; CSEN Wound needle - 16G x 80 mm;

Wound needle filter (blue); Wound catheter filter (white). Packed, Sterilized,

CE marked.

CSEN 194.1: - **ELDOR Wound Catheter kit** - OD 1 mm x L 1000 mm, 20

lateral holes, 10 mm apart, closed tip; Wound catheter adapter (Tuohy Borst);

Wound catheter thread assist guide; CSEN Wound needle - 16G x 80 mm;

Wound needle filter (blue); Wound catheter filter (white).). Packed, Sterilized,

CE marked.

CSEN 194.2: - **ELDOR Wound Catheter kit** - OD 1 mm x L 1000 mm, 10

lateral holes, 10 mm apart, closed tip; Wound catheter adapter (Tuohy Borst);

Wound catheter thread assist guide; CSEN Wound needle - 16G x 200 mm;

Wound needle filter (blue); Wound catheter filter (white). Packed, Sterilized,

CE marked.

CSEN 194.2D: - **ELDOR Wound Catheter kit** - OD 1 mm x L 900 mm, 10

lateral holes, 10 mm apart, closed tip; Wound catheter adapter (Tuohy Borst); Wound catheter thread assist guide; CSEN Wound needle - 16G x 200 mm; Wound needle filter (blue); Wound catheter filter (white). Packed, Sterilized, CE marked.

CSEN 194.2 - Cat: - ELDOR Wound Catheter kit - OD 1 mm x L 1000 mm, 10 lateral holes, 10 mm apart, closed tip; In Bulk.

CSEN 194.3: - ELDOR Wound Catheter kit - OD 1 mm x L 1000 mm, 10 lateral holes, 10 mm apart, closed tip; Wound catheter adapter (Tuohy Borst); Wound catheter thread assist guide; Wound needle filter (blue); Wound catheter filter (white). Packed, Sterilized, CE marked.

CSEN 194.4: - CSEN Wound needle - 16G x 200 mm. Packed, Sterilized, CE marked.

CSEN 194.5: - ELDOR Wound Catheter - OD 1 mm x L 1000 mm, 20 lateral holes, 10 mm apart, closed tip. Packed, Sterilized, CE marked.

CSEN 194.6: - ELDOR Wound Catheter kit - OD 1 mm x L 1000 mm, 20 lateral holes, 5 mm apart, closed tip; Wound catheter adapter (Tuohy Borst); Wound catheter thread assist guide; CSEN Wound needle - 16G x 200 mm; Wound needle filter (blue); Wound catheter filter (white). Packed, Sterilized, CE marked.

CSEN 194.7: - ELDOR Wound Catheter kit - OD 1 mm x L 1000 mm, 40 lateral holes, 5 mm apart, closed tip; Wound catheter adapter (Tuohy Borst); Wound catheter thread assist guide; CSEN Wound needle - 16G x 200 mm; Wound needle filter (blue); Wound catheter filter (white). Packed, Sterilized, CE marked.

